

iMID 2015 Display Special Exhibition

Materials, Components and Instrument

Company Name	Surface Electro Optics(S.E.O)
Address	112,Sanup-ro Kwonsungu Suwon-City Kyuggido 441-813, Korea
President	Park Jang ho
Tel	82.31.298.9561~4
FAX	82.31.298.9565
Email	seo@s-eo.com
Homepage	www.s-eo.com
Exhibitor Introduction (Within 200 words)	<p>SEO opened a new stage of Contact Angle Analysis Technology as a result of a lot of effort for the localization of Surface Chemistry Analysis instruments. In addition, ever since the SEO Phoenix Series came into the market, we have been enjoying a very good reputation for quality and service. Since SEO is based on the development of technology of our own, the instruments we manufacture have more competitive price than those made by international companies and compare very favorably in quality</p>
Exhibit Description (Within 200 words)	<p>Contact Angle Analyzer is the equipment for measuring the changes in the properties of a solid surface, A device for evaluating the surface relatively easily and accurately compare change. Liquid according to the method of surface evaluation The surface energy or adhesive strength of each adhesive or solid wetting spreadability of the solid, and the like can be evaluatedA computer controlled contact angle analyzer makes it easy to measure static and dynamic contact angle, along with surface tension and surface energy.</p>
Exhibit Product	<p>Contact Angle Analyzer (Model: Phoenix 300 Touch)</p>
Company logo	